(EN)

**Instructions for use**

Bleaching gels developed solely for in-office bleaching in the dental practice. A double-cartridge for instant mixing (4:1) contains bleaching gel and activator. The mixed gel contains 6% hydrogen peroxide. May only be sold to dentists. The first round of treatment must be performed by a dentist, or under direct supervision of a dentist, if equivalent safety can be guaranteed. May thereafter be handed over for continued treatment. Remove the cap from the cartridge and replace with the mixing tip. Press the gel out of the cartridge until a uniformly colored gel emerges from the mixing tip. Discard the inhomogeneous material which emerged initially. Storage at higher temperatures for a longer period of time can cause the development of isolated bubbles within the material. This may lead to inhomogeneous, liquid mixtures, especially when there is less than a quarter of the material left inside the double cartridge. In this case, do not use the remaining material.

**The ready-mixed bleaching gel irritates the skin and the mucous membranes, and it causes severe eye damage. Always wear protective gloves and goggles when using the material! Avoid contact of the material with eyes, skin and mucous membranes.**

**Indications**

- Brightening of devital teeth
- Brightening of pathologically discolored teeth as an alternative to crown or veneer
- Treatment of discolorations caused by pharmaceuticals (e.g. Tetracycline)
- Aesthetic reasons

**Contraindications**

- Exposed dentine
- Defective filling margins
- Patients with hypersensitivities
- Allergies or incompatibilities with the ingredients of the product
- Pregnancy or breast-feeding (based upon a lack of experience)

Many patients may experience slight sensitisation in their teeth during or shortly after treatment, but this normally subsides after a short time.

**Stop treatment immediately if undue pain is being experienced.**

**Preparation**

Determine the shade of the tooth by comparison with a Vita® shade guide (photo 1). For assessment, re-arrange the Vita® shade guide as listed below.

Values arranged in accordance with the shades:

A1 A2 A3 A3.5 B1 B2 B3 B4 C1 C2 C3 C4 D1 D2 D3 D4

Values arranged in accordance with the brightness:

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B2 to B1 = increase of value by 2

A3 to A2 = increase of value by 4

**Step by step instructions:**

1. The patient cleans the teeth with Dental Cleanser, a cleaning cloth to put on the finger which removes saliva and coatings from the teeth. The patient should not close the mouth before inserting the tray, as saliva counteracts the bleaching process.
2. Apply a thin layer of gel on the foam strip of the tray and spread it into an even layer.
3. Apply vaseline to the patient's lips and then place the tray into the patient's mouth. The foam strip should be against the teeth and the lips should rest outside of the tray. Be sure to inform that the patient can and should swallow their saliva during the treatment.
4. Place the lamp in front of the mouth. The custom made tray has a magnifying glass in front, spreading the light in the mouth. The teeth are illuminated for 20 minutes at a time.
5. After 20 minutes, take out the tray, wipe dry with tissue and apply a new layer of gel. The patient can rinse the mouth during this time.
6. Dry the teeth with tissue, then again place the tray in the mouth and illuminate for 20 minutes. If necessary, this process is repeated once more.
7. When the light treatment is complete, remove the disposable tray. For sanitary reasons, the gel tip is removed and disposed after each patient.

**Note:** Composite fillings or ceramic and metal restorations cannot be

bleached. Such fillings or restorations should be replaced about two weeks after the treatment.

**Precautionary instructions**

- ! **The patient and the practitioner should wear safety goggles when light treatment is in progress.**
- ! **Do not use an excess amount of gel, as this may result in excess gel ending up on the gums causing irritation.**
- ! **Remove any excess gel on the lips or on the mucosal tissue in the mouth.**

**After the treatment**

A sealing of the teeth with NanoSeal Total+ can be offered to patients who wish to be able to eat and drink anything right after treatment. The solution contains nanostructured fluoride which rebuilds the enamel immediately. NanoSeal Total+ is applied to dry and clean teeth after the teeth whitening treatment is completed.

- ! Without sealing, the following instructions are applicable: 72 hours after completion of the treatment, the patient should avoid acidic foods and drinks (cola, citrus drinks and fruit), as these can cause hypersensitivity. Also highly colored drinks (coffee, tea, red wine) and food (chocolate, tomatoes, curry powder, oranges, etc.) and tobacco should be avoided during the same time period, as this may compromise the results.

**Shelf-life/Storage**

Gel should not be used after its expiration date, as a sufficient bleaching performance can then no longer be guaranteed. **Store the bleaching gel in a refrigerator (3-9°C), but do not freeze.**

**Product contains hydrogen peroxide, potassium nitrate, sodium fluoride, sodium hydroxide, glycol derivatives, thickening agent and dye.**

**Possible sources of error**

Bleaching result	Repeat bleaching process (up to five treatment cycles in one sitting). Additional bleaching process can be performed after one week.
Patient experiences excessive discomfort or pain during treatment	Gel is in contact with exposed dentine or an untight filling (check before treatment) Patient has very sensitive teeth (depends on the individual tooth structure)
White stains on the enamel after bleaching	Local demineralization or dehydration: this is reversible and disappears after fluoridation
Gingiva appears white after bleaching	This is caused by an insufficient or untight gingiva protection, but is also reversible.
Gel protruding from the cartridge is not uniformly colored	Caused by air bubbles in the gel, which appear if the gel is stored at higher temperatures than recommended or if it is stored overlong peroxide residue
Composite enamel adhesion decreased after treatment	Wait 1-2 weeks before replacing restorations

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Bestem farge på tannen ved å sammenligne med en Vitapan© fargeskala (bilde 1). For korrekt vurdering, omarbeid Vita© fargeskala som listet opp nedenfor.