# VICRYL™ Polyglactin 910 Suture

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# VICRYL™

## Polyglactin 910 Suture

### Sterile Synthetic Absorbable Surgical Suture

### Description

VICRYL™ Suture is a sterile, coated/uncoated, braided/monofilament, synthetic, absorbable, surgical suture composed of ≥ 87.8 weight percent Polyglactin 910, a copolymer made of 90 weight percent glycolide and 10 weight percent L-lactide. The empirical formula of the copolymer is  $(C_2H_2O_2)_m(C_3H_4O_2)_n$ .

VICRYL™ Suture is available as a coated, braided suture; some sizes are also available as uncoated braided suture (size 9-0) and as an uncoated monofilament suture (sizes 9-0 and 10-0). VICRYL™ Suture has a coating with up to 12.0 weight percent of a mixture composed of approximately equal parts of calcium stearate and Polyglactin 370, a copolymer composed of 30 weight percent glycolide and 70 weight percent lactide. Polyglactin 910 copolymer and Polyglactin 370 with calcium stearate have been found to be nonantigenic, nonpyrogenic and elicits only a slight tissue reaction during absorption.

VICRYL™ Suture is available undyed and dyed with up to 0.20 weight percent D&C Violet No. 2 (Color Index 60725) to enhance visibility in the surgical field.

VICRYL™ Suture is available in a range of gauge sizes and lengths, non-needled or attached to needles of various types and sizes, and in presentations as described in the How Supplied section.

VICRYL™ Suture complies with the requirements of the European Pharmacopoeia (Ph. Eur.) for Sterile Synthetic Absorbable Braided Sutures and United States Pharmacopoeia (USP) for Absorbable Surgical Sutures, except for a slight oversize in diameter for the gauges as shown in the table below.

USP Suture Size	Metric Suture Size	Maximum Oversize
9-0 (uncoated braid only)	0.3	0.010 mm
8-0	0.4	0.020 mm
7-0	0.5	0.030 mm
6-0	0.7	0.050 mm
5-0	1	0.050 mm
4-0	1.5	0.050 mm
3-0	2	0.090 mm
2-0	3	0.060 mm
0	3.5	0.100 mm
1	4	0.100 mm
2	5	0.100 mm
3	6	0.100 mm

The European Pharmacopoeia recognizes units of measure Metric and Ph. Eur. Sizes as equivalent, which is reflected on the labeling.

VICRYL™ Suture is intended for use only by healthcare professionals who are trained in surgical suturing techniques.

Clinical benefits expected from soft tissue approximation and/or ligation are to promote wound healing by primary intention and to avoid bleeding from vessels during and after surgical procedure.

A summary of safety and clinical performance can be found at the following link (upon activation):

<https://ec.europa.eu/tools/eudamed>.

## Indications / Intended Use

VICRYL™ Suture is indicated for use in general soft tissue approximation and/or ligation, including use in ophthalmic surgery, but not for use in cardiovascular and neurological tissues.

## Contraindications

This suture, being absorbable, should not be used where extended approximation of tissue under stress is required.

## Warnings

The safety and effectiveness of VICRYL™ Suture has not been established in the following areas: cardiovascular and neurological tissues.

Healthcare professionals should be familiar with surgical procedures and techniques involving absorbable sutures before employing VICRYL™ Suture for wound closure, as risk of tissue separation / wound dehiscence leading to impaired healing may vary with the site of application and the suture material used. Healthcare professionals should consider the in vivo performance (under the Performance / Actions section) when selecting a suture.

Acceptable surgical practice should be followed for the management of contaminated or infected wounds.

As this is an absorbable suture material, the use of supplemental non-absorbable sutures should be considered by the healthcare professional in the closure of the sites that may undergo expansion, stretching or distension, or that may require additional support.

As an absorbable suture, this suture may act transiently as a foreign body. As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation. Like all foreign bodies, VICRYL™ Suture may potentiate infection.

Do not resterilize/reuse. Reuse of this device (or portions of this device) may create a risk of product degradation that may result in device failure and/or cross-contamination that may lead to infection or transmission of blood-borne pathogens to patients and healthcare professionals.

Discard unintentionally opened / partially used / used devices and packages.

## Precautions

As with any suture material, adequate knot security requires the standard surgical technique of flat, square ties with additional throws as warranted by surgical circumstance and the experience of the healthcare professional. Improper knotting or damage to

suture during use can cause extended surgery time or treatment failure and additional surgery. The use of additional throws may be particularly appropriate when knotting monofilament sutures.

Skin sutures that must remain in place longer than 7 days may cause localized irritation and should be snipped off or removed as indicated by the healthcare professional. Subcuticular sutures should be placed as deeply as possible to minimize the erythema and induration normally associated with the absorption process.

Under some circumstances, notably orthopedic procedures, immobilization of joints by external support may be employed at the discretion of the healthcare professional.

Consideration should be taken in the use of absorbable sutures in tissues with poor blood supply as suture extrusion and delayed absorption may occur.

This suture may be inappropriate in elderly, malnourished or debilitated patients, or in patients suffering from conditions that may delay wound healing.

In handling this or any other suture material, care should be taken to avoid damage. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

Care should be taken to avoid damage when handling surgical needles. Grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the attachment end to the point. Grasping in the point area could impair the penetration performance and cause fracture of the needle. Grasping at the attachment end could cause bending or breakage. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking.

Healthcare professionals should exercise caution when handling surgical needles to avoid inadvertent needle stick injury that may result in transmission of blood-borne pathogens from contaminated needles. Broken needles may result in extended or additional surgeries or residual foreign bodies. In the event of a product malfunction before use, such as a bent, broken or detached needle, or suture damage, the product should be discarded and a new one obtained to begin the procedure. In the event of product malfunction during use, it is up to the discretion of the healthcare professional whether to continue or discontinue usage of the product and how to complete the procedure.

## Adverse Reactions / Undesirable Side Effects

Adverse reactions associated with the use of this device include tissue separation / wound dehiscence leading to impaired healing or failure to provide adequate wound support in the closure of the sites where expansion, stretching, or distension occur, and in elderly, malnourished or debilitated patients or in patients

suffering from other conditions that may delay wound healing. Other associated adverse events include calculus formation in urinary or biliary tracts when prolonged contact with salt solutions such as urine or bile occurs, transient inflammatory tissue reaction, and delayed absorption in tissue with poor blood supply. Transient local irritation and suture extrusion may occur at the wound site when skin sutures are left in place for greater than 7 days. Like all foreign bodies, VICRYL™ Suture may potentiate infection.

Suture breakage may result in hemorrhage. Suture material cutting or tearing through the tissue at the time of suture tensioning may cause soft tissue injury. Improper knotting or damage to suture during use can cause extended surgery time or treatment failure and additional surgery.

Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of blood-borne pathogens.

Healthcare professionals should convey adverse reactions, undesirable side effects and risks associated with the product and the procedure to the patient and advise the patient to contact a healthcare professional in case of any deviation from the normal postoperative course.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the country-competent authority.

### **Magnetic Resonance Imaging (MRI) / Carcinogenic, Mutagenic, Toxic to Reproduction (CMR) / Endocrine Disrupting (ED) Safety Information**

VICRYL™ Suture is MR safe. No known CMR Category 1a/1b and ED substances are present at >0.1%. Category 1a/1b are defined as known or presumed human carcinogen (H340), mutagen (H350) or reproductive toxicant (H360), based on human evidence and animal studies.

### **Application / Instructions for Use**

Sutures should be selected and implanted depending on patient condition, surgical experience, surgical technique and wound characteristics. Discard needles in “sharps” containers. Discard unintentionally opened / partially used / used devices and packages according to your facility’s policies and procedures concerning biohazardous materials and waste.

### **Performance / Actions**

VICRYL™ Suture elicits an initial, minimal inflammatory reaction in tissue and ingrowth of fibrous connective tissue. Progressive loss of tensile strength and eventual absorption occurs by means

of hydrolysis, where the copolymer degrades to glycolic and lactic acids, which are subsequently absorbed and metabolized in the body. Absorption begins as a loss of tensile strength followed by a loss of mass. All of the original tensile strength is lost by five weeks post implantation. Absorption is essentially complete between 56 and 70 days. Implantation studies in rats indicate that VICRYL™ Suture retains the percentage of original tensile strength as follows:

<b>Days Implantation</b>	<b>Approximate Original Strength Remaining</b>
14 Days	75%
21 Days USP 6-0 (Metric 0.7) and larger	50%
21 Days USP 7-0 (Metric 0.5) and smaller	40%
28 Days USP 6-0 (Metric 0.7) and larger	25%

### **Sterility**

VICRYL™ Suture is sterilized by ethylene oxide gas. Do not resterilize. Do not use if package is opened or damaged.

### **Storage**

No special storage conditions required. Do not use after expiry date.

### **How Supplied**

Please note that not all sizes are available in all markets. Please contact your local sales representative for size availability.

VICRYL™ Suture is available as sterile strands in sizes 10-0 through 3 and 6 (metric sizes 0.2-6.0 and 8.0), in a variety of lengths, with and without permanently attached needles.

VICRYL™ Suture is also available as a monofilament in sizes 9-0 and 10-0 (metric sizes 0.3 and 0.2).

VICRYL™ Suture is also available as sterile strands attached to CONTROL RELEASE™ Needles that enable the needles to be pulled off instead of being cut off.

VICRYL™ Suture is available in one, two or three dozen units per box.

### **Traceability**

The following specific information can be found on the device packaging label: Catalogue number, Batch Code, expiry and manufacturing date, manufacturer name, address and the website and a Unique Device Identification bar code with the Global Trade Item Number information.

## Symbols Used on Labeling

	Catalogue number		Dyed - Absorbable - Braided - Uncoated - Suture
	Medical Device		Undyed - Absorbable - Braided - Coated - Suture
	Do not use if package is damaged and consult instructions for use		Dyed - Absorbable - Braided - Coated - Suture
	Caution		Dyed - Absorbable - Monofilament - Uncoated - Suture
	Do not re-use		CONTROL RELEASE™
	Do not re-sterilize		Multi-strand
	Single sterile barrier system with protective packaging inside		
	Sterilized using ethylene oxide		
	Batch code		
	Date of manufacture		
	Use-by date		
	Unique Device Identifier		
	Manufacturer		
	Authorized Representative in the European Community/European Union		
	Packaging unit		
	PEEL/LIFT/OPEN HERE		
	TEAR HERE		
 <small>www.e-ifu.com EU: +800 8888 2020 EU: +32 2 4037222</small>	Consult instructions for use or consult electronic instructions for use. EU: Call paper on demand help desk to get paper copies free of charge within 7 days.		
	MR Safe		

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AB: Basılı kopyaları 7 gün içinde ücretsiz olarak talep üzerine almak için danışma hattını arayın.
Johnson & Johnson Medical GmbH
Ethicon, LLC