O-Bite

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Vinylpolysiloxane, Hydrogenpolysiloxane, Füllstoffe, Pigmente, Additive, Platinkatalysator

Handelsformen

Packung	
2 Kartuschen à 50 ml Paste, 12 Automix-Tips	REF 909765*
Vorteilspackung	
8 Kartuschen à 50 ml Paste, 50 Automix-Tips	REF 909767*
Zubehör	
50 Automix-Tips	REF 909202*
50 Conturation-Tips	REF 909203*

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Packung	
4 Kartuschen à 50 ml Paste, 15 Automix-Tips	REF 999765
Zubehör	
50 Automix-Tips	REF 909202

Instructions for Use

English

Product description

Bite registration material, A-silicon-based, with high end hardness.

Indications

- Bite registration
- Material for application on bite forks for use with facebows
- Key material, e.g. for intraoral incisal pin records

Contraindications

Do not use the material in cases of allergies to any of the ingredients or, as can occur in rare cases, in the event of contact allergies.

Basic safety instructions

For dental use only!
Keep out of reach of children!

Interactions

Never allow the material to come into contact with organic solvents or liquids containing solvents as this may result in a swelling of the material and may therefore lead to inaccuracies.

Avoid contact with retraction solutions, because in some instances residues of retraction solutions can increase the setting time. Remove retraction solutions by spraying and rinsing thoroughly.

Do not mix the material with condensation curing silicone or polyether impression materials.

Avoid contact between the impression material and latex gloves because these prevent the setting process. In the event of prior contact between the hands and latex gloves, wash hands thoroughly. It is recommended that PVC or PE based gloves be used.

Using the Automix cartridge

Insert cartridge

1. Push up the lever [Fig. 1, Pos. A] at the back of the Automix dispenser and pull back the slide [Fig. 1, Pos. B] completely.
2. Lift up the plastic lever [Fig. 1, Pos. C], insert the cartridge and lock it using the plastic lever [Fig. 1, Pos. C].

Please note: Make sure that the notches on the cartridge line up with those on the Automix dispenser.

Attach mixing tip

1. Remove the cartridge cap [Fig. 1, Pos. D] or used mixing tip by turning it 90° counterclockwise, and discard it.
2. Attach a new mixing tip [Fig. 1, Pos. E].

Please note: In order to attain an optimum mixing result, DMG recommends the use of the mixing tips available from DMG. All mixing tips are also available as refill packs.

Please note: Make sure that the notches on the mixing tip line up with those on the cartridge.

3. Turn the mixing tip clockwise 90° to lock it.

Applying the material

When using a cartridge for the first time, squeeze out a pea-sized amount of material and discard it.

The material is mixed in the mixing tip by pulling the trigger [Fig. 1, Pos. F] of the Automix dispenser and it is then applied directly.

Please note: After the application, leave the used mixing tip on the cartridge to cap unused material!

Remove cartridge

1. Push up the lever [Fig. 1, Pos. A] at the back of the Automix dispenser and pull back the slide [Fig. 1, Pos. B] completely.
2. Lift up the plastic lever [Fig. 1, Pos. C] and remove the cartridge.

Recommended use

Create bite registration

1. Apply the bite registration material directly to the occlusal surfaces. If necessary, apply to prepared areas first.

Please note: It is usually sufficient to apply the material to the teeth of the upper or lower dental arch. In the area of the front teeth, it is recommended that the material be applied to both dental arches. A material layer approx. 5 mm thick should be sufficient.

2. Allow patient to adopt the desired occlusion.
3. After 1:30 min., carefully remove the finished bite registration from the mouth.

Process bite registration

Please note: The material has a high final hardness and can be processed with a rotating instrument (e.g. grinder). Immediately after removal, excess can also be removed with a scalpel.

1. Remove excess material extending beyond the cast chewing surface areas.
2. Remove shreds of material which extend into the approximal areas.

Disinfection of the bite registration

Disinfect bite registration with standard disinfection solutions.

Storage of the bite registration

Special storage conditions are not required. The bite registration can be stored for a long period of time.

Disinfection of the Automix dispenser

The Automix dispenser is reusable after disinfection, although it should be replaced if visible signs of wear appear.

Disinfection can be carried out using a commercially available dipping disinfectant. DMG recommends using only disinfecting solutions listed by the RKI (Robert Koch Institute). The Automix dispenser is also autoclavable.

Technical data

Total working time*	≤ 0:30 min
Recommended time in mouth	≥ 1:30 min
Shore-A hardness (1 h)	≥ 93

*The times apply at a room temperature of 23 °C and a normal relative air humidity of 50%. Higher temperatures shorten, lower temperatures increase this time.

Storage and shelf life

Store at room temperature (15 – 25 °C/59 – 77 °F) in a dry place.
Do not store in the fridge or freezer.
Protect against exposure to direct sunlight.
After the application, leave the used tip on the cartridge to cap unused material.
Do not use after the expiration date.

Composition

Vinyl polysiloxanes, hydrogen polysiloxanes, fillers, pigments, additives, platinum catalyst

Available product sizes

Packaging	
2 cartridges @ 50 ml paste, 12 Automix-Tips	REF 909765*
Economy pack	
8 cartridges @ 50 ml paste, 50 Automix-Tips	REF 909767*
Accessories	
50 Automix-Tips	REF 909202*
50 conturation tips	REF 909203*

*Not available in the United States.

Only available in the United States

Packaging	
4 cartridges @ 50 ml paste, 15 Automix-Tips	REF 999765

Accessories	
50 Automix-Tips	REF 909202

Federal law (USA) restricts this device to sale by, or on the order of a dentist, or other practitioner licensed by law of the state in which he or she practices to use or order the use of this device.

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