

# UNIVERSAL PROSTHETIC PRIMER FOR MULTIPLE CEMENT AND REPAIR INDICATIONS

CLEARFIL™ CERAMIC PRIMER PLUS is a universal prosthetic primer for all cement and repair indications. The product incorporates both Kuraray Noritake's original MDP monomer and a silane monomer. This allows for excellent bond strength to all silica-based ceramics, zirconia, composites and even metals.



CLEARFIL™ CERAMIC PRIMER PLUS is a universal, apply, dry and proceed primer for use with most restorative materials. No waiting time after applying! In adhesive cementing procedures it is absolutely 100% compatible with any resin-based cement, regardless the brand. For example; RelyX (Solventum\*), Calibra (Dentsply Sirona\*) G-CEM (GC\*), Variolink (Ivoclar Vivadent\*), BiFix (VOCO\*), Choice (Bisco\*), Maxcem (Kerr\*).

## CLEARFIL™ CERAMIC PRIMER PLUS is indicated for:



Surface treatment of prosthetic restorations made of ceramics\*\*, hybrid ceramics, composite resin or metal.



Intraoral repair of fractured restorations made of ceramics\*\*, hybrid ceramics, composite resin or metal.

#3637-EU CLEARFIL™ CERAMIC PRIMER PLUS 4ml

\* No registered trademarks or trademarks of KURARAY CO., LTD.\*\* Metal Oxide, glass ceramics.

