

CALASEPT® Plus

Calcium Hydroxide

”The most effective calcium hydroxide thanks to high concentration and optimal calcium release”



NEW NEW NEW

CALASEPT® Plus

- Calcium Hydroxide

CALASEPT Plus has a very high content of calcium hydroxide (> 41%) giving a high concentration of hydroxyl ions. These generate the high pH value of 12.4 which has a pronounced bactericidal effect, i.e. killing of bacteria. CALASEPT Plus is radio-opaque.

High concentration of calcium hydroxide means long-lasting because the paste can release calcium ions for a long time. CALASEPT Plus is a ready-to-use paste in air tight syringes for direct application through the Flexi-tip.

■ FACT	● BENEFIT
High concentration, >41%	Powerful, long-lasting effect
Pure calcium hydroxide	Optimal ion release
Smooth consistency + Flexi-tip	Precise and deep access
Pre filled syringe	Direct, quick and easy application
Multi purpose	Five indications – only one product

COMPOSITION
Calcium hydroxide 41 g
Barium sulphate 5 g
Sterile isotonic saline solution pH 12.4



Flexi-tip (autoclavable) in combination with the smooth consistency allows precise and deep application.




PACKAGING:

CALASEPT Plus 4U
4 syringes each containing 1,5 ml CALASEPT Plus
20 CALASEPT Plus Flexi-tips, luer lock

CALASEPT Plus 1U
1 syringe containing 1,5 ml CALASEPT Plus
5 CALASEPT Plus Flexi-tips, luer lock

CALASEPT PLUS Flexi-tip
100 CALASEPT Plus Flexi-tips, luer lock



Box 1082, SE-262 21 Ängelholm, Sweden
Int tel: (+46) 431 44 33 60 - Fax: (+46) 431 44 33 99
mail@nordiskadental.se - www.nordiskadental.se
Certified acc. to ISO 9001 and ISO 13485