


Dual Cure
Calibra[®]
Esthetic Resin Cement

Finish Strong.



Easy
clean up!

NEW!



Calibra[®]
Esthetic Resin Cement

For better dentistry

DENTSPLY

Finish Strong.

Calibra[®] Resin Cement is available in a dual cure Automix Syringe. You can now take advantage of a complete shade selection with shade stable[™] chemistry in the delivery system of your choice.

✓ Easy clean up

Light cure excess cement for 10 seconds to initiate a gel phase that makes clean up fast and simple.



✓ Full spectrum of esthetic shading possibilities



Translucent



Light



Medium



Dark



Opaque

Displayed colors may not perfectly represent the actual product

DENTSPLY's Shade Stable[™] chemistry gives you the confidence that the shades will not shift over time.

Esthetic Resin Cement



For virtually any indirect restoration

- Ceramic, porcelain, composite inlays/onlays, Maryland bridges, veneers and crowns
- All metal crowns, bridges, inlays/onlays including precious, semi-precious and non-precious metals
- PFM (porcelain fused to metal) crowns and bridges
- Cementation of prefabricated and cast posts



Multi-Unit Bridge



Cercon® Zirconia



All-Ceramic Crowns & Bridges




Captex® PFM



Long Lasting Strength


Shear Bond Strength when adhesive is not light cured¹

Shear bond strength [MPa] after 24 hours at 37°C



Marginal Integrity of ceramic inlays after thermal cycling³

Marginal Integrity [%]



Results:

XP BOND® + Self Cure Activator + Calibra® is equal to, or better than other leading adhesive cementation systems.

¹ Internal Data

² Not a registered trademark of DENTSPLY International, Inc.

³ Source: Frankenberger R, et al., Luting of ceramic inlays in vitro: Marginal quality of self-etch and etch-and-rinse adhesives versus self-etch cements, Dental Mater (2007), doi: 10.1016/j.dental.2007.04.003

Calibra®
Dual Cure

Ordering Information

Order-No.

Calibra® Automix

1 Automix Syringe (4.5 g)



- 607084 ○ Translucent
- 607081 ○ Light
- 607082 ○ Medium
- 607083 ○ Dark
- 607085 ○ Opaque

Mixing Tips Refill

50 Standard Mixing Tips



607086

Calibra® Base Refill Pack

1 Syringe (2 g)



- 607064 ○ Translucent
- 607061 ○ Light
- 607062 ○ Medium
- 607063 ○ Dark
- 607065 ○ Opaque

Calibra® Catalyst Refill Pack

1 Syringe (2 g)



- 607068 Regular Viscosity
- 607069 High Viscosity

Operatory Kit

2 Syringes Base: 1x Translucent, 1x Light (2 g each),
 2 Syringes Catalyst: 1x Regular Viscosity, 1x High Viscosity (2 g each),
 1 Bottle Prime&Bond® NT (4.5 ml), 1 Bottle Self-Cure Activator (4.5 ml),
 1 Syringe Calibra® Silane Coupling Agent (3 ml), Dental Brush Handle,
 Brush Tips, 3 Well Mixing Trays, Mix Stick, Silane Tips,
 1 Illustrated Technique Guide



607059

2-Part Complete Kit

Kit I: 5 Syringes Calibra® Base, 2 g each: Light, Medium, Dark, Translucent, Opaque
 2 Syringes Calibra® Catalyst, 2 g each: Regular and High Viscosity
 1 Bottle Prime&Bond® NT (4.5 ml), 1 Bottle Self-Cure Activator (4.5 ml),
 1 Syringe Tooth Conditioner Gel (3 ml), Dental Brush Handle, Brush Tips,
 3 Well Mixing Trays, Mix Stick, Application Needles, Illustrated Technique Guide.
 Kit II: 5 Syringes Calibra® Try-In Paste (1.8 g each): Light, Medium, Dark, Translucent, Opaque
 1 Syringe Calibra® Silane Coupling Agent (3 ml), Amber Cover, Mixing Pad, Silane Tips.



607060

Try-In Accessory Pack

5 Syringes Try-In Paste – 1 each shade: Translucent, Light, Medium, Dark, Opaque

607058

DENTSPLY DeTREY GmbH

De-Trey-Str. 1
 78467 Konstanz
 GERMANY
 www.dentsply.eu
 Phone +49 (0) 75 31 5 83-0

DENTSPLY Ltd

Building 1
 Aviator Park, Station Road
 ADDLESTONE, KT15 2PG
 UNITED KINGDOM
 Phone +44 (0) 19 32 85 34 22

For better dentistry

DENTSPLY